

## QC-603 Manual Specimen Punching Machine

### Standard

---

### Description

---

This machine can be with variety of cutter to cut the specimen or standard shape of specimen for do tensile strength, tear, elongation, aging and sample...etc. of test. It is suitable cutting material: paper, rubber, soft plastic and foam...etc of soft materials.

### Industry

---

Rubber, Plastic

### Appearance

---



## Specification

---

|                      |   |
|----------------------|---|
| Model                | QC-603  |
| Suit for material    | rubber, soft plastic and leather, etc.                  |
| Pressing Board area  | 70×150mm  |
| Accessories          | Cushion board one piece                                 |
| Operation method     | Steering wheel to descend cutter and punch specimen     |
| Optional accessories | Cutter, please provide specimen dimension to customized |
| Dimension            | 42 × 42 × 75 cm(W×D×H)                                  |
| Weight               | 65 kg   |

## Cutter type (Optional)

---



A type



B type

## Related product

---

QC-601T Manual Thermos Press Forming Machine

QC-601A Manual Thermos Press Forming Machine

QC-601B Manual Thermos Press Forming Machine

QC-602A Hydraulic Hot Forming Press Machine

QC-602B Hydraulic Press Forming Machine

QC-603A & QC-603C Pneumatic Specimen Punch

QC-603G Hardness Specimen Cutting Device

\*Cometech reserves ultimate modification right of product specification.

---



COMETECH TESTING MACHINES CO.,LTD.

No.49, Ln, 39, Liancun Rd, Fengyuan Dist, Taichung City 420, Taiwan

Tel/886-4-2536881 Fax/886-4-2536883

e-mail/sales@come-tech.com.tw [www.come-tech.com.tw](http://www.come-tech.com.tw)